

**Easily Accommodates Patients with Larger Calves.  
800 lbs (363 kg) Patient Weight Capacity.**

*Note: Stirrups are Sold in Pairs*

The Ultrafins Stirrup is an upgrade of the Yellofins® Stirrup. The reinforced rod and joint assembly and the flexible boot make the system capable of positioning a patient weighing up to 800 lbs (363 kg). Though designed for large patients, the Ultrafins Stirrup is still appropriate for use with the general patient population.

- **Patient Weight Capacity: 800 lbs (363 kg):** Support lithotomy positioning for very large patients.
- **Flexible Boot:** Boot comfortably conforms to the calf shape of the patient while providing support.
  - Fin design helps reduce pressure under the popliteal fossa and superficial peroneal nerve.
  - Flexible boot allows these stirrups to be standardized throughout the facility.
- **Floating Boot:** Self adjusting boot minimizes pressure on the calf, when moving the stirrup.
- **Lift-Assist™ Technology:** Allows easy movement of the leg when placing it in the desired position.
- **Ultrafins Clamshell Pad:** Completely encapsulates the foot, ankle and calf for comfortable and secure positioning for both large and small patients.
- **Squeeze Grip Handle:** Provides easy intraoperative adjustment without compromising the sterile field. Simply release the handle to secure the leg holder in all directions.
- **Lithotomy Range:** Set the stirrup in any position between +84° to -33° lithotomy.
- **Abduction Range:** Set the stirrup in any position between +25° to -9° abduction.
- **Lithotomy & Length Indicators:** Help ensure precise positioning.

#O-UFAS Ultrafins Stirrups

<b>What's Included</b>	 <p>Ultrafins Stirrups #O-UFAS 1 Pair</p>	 <p>Ultrafins Clamshell Pad #O-AMA65S0.6S 1 Pair</p>	
<b>What's Required</b>	 <p>Easy Lock Blade Clamp #A-40040 (US), #A-40041 (EU), #A-40042 (UK), #A-40043 (Jpn), #A-40044 (Den), 2 Required, Specify Rail Type</p>	 <p>Safety Drape #A-22000 Case of 10 Pairs</p>	
<b>What's Optional</b>	 <p>Easy Catcher Disposable System (Case of 10) (Sterile) (Non-Sterile) (See Page 85)</p>	 <p>Stirrup Cart #A-30015</p>	 <p>Extended Warranty #A-XWARUF Available in US Only</p>



**Ultrafins® Stirrups**  
Supports up to 800 lbs (363 kg)



**Patient Weight Capacity:**  
800 lbs (363 kg)



**Clamshell Pad Shown**