

**Easy to Use. Fully Disposable Single Use System**

*The All-in-One Fluid Containment System*

The Allen® Easy Catcher is designed to catch and drain fluid during a surgical procedure. It attaches to the side rails of an OR table via two integral rail clamps. The Easy Catcher is a single-piece system that contains a fluid drain drape with an integral flexible support band. The band supports the drape vertically for optimal drainage. The support band also flexes to accommodate the doctor’s mid-section during procedures.

The all-in-one system is fully disposable making it more convenient and quicker to setup and dispose of. In addition, the Easy Catcher has a removable tube adaptor that enables it to connect to the Allen Uro Collection Container and a variety of fluid waste management systems.

- Saves clean-up time
- Fully disposable
- Flexible integrated band
- Easy set-up and removal
- Integral clamps for easy set-up
- Single-piece system includes: buttocks flap, mesh trap, collapsible hose and removable tube adaptor
- Maximize infection control with disposable system
- Attaches to the rails of any OR table
- Eliminates the need to clean system band for repeated use
- No more lost or damaged bands
- Removable tube adaptor connects to a variety of fluid waste management systems

**Ideal for use with a variety of urology and gynecology procedures including:**

- Cystoscopy
- Ureteroscopy
- Urologic Oncology
- Pediatric Urology Procedures
- Calculi (urinary tract stones) Procedures
- Neurourology Procedures

#O-EC51 Easy Catcher Disposable System: STERILE (Case of 10) (US)

#O-EC52 Easy Catcher Disposable System: NON-STERILE (Case of 10) (US)

#O-UCC20S Allen Uro Collection Container, 20L with Solidifier (Case of 10)

**Specify Rail Type**

- |                            |                                |
|----------------------------|--------------------------------|
| #O-EC51 (US) Sterile       | #O-EC52 (US) Non-Sterile       |
| #O-EC51-EU (EU) Sterile    | #O-EC52-EU (EU) Non-Sterile    |
| #O-EC51-UK (UK) Sterile    | #O-EC52-UK (UK) Non-Sterile    |
| #O-EC51-JP (Japan) Sterile | #O-EC52-JP (Japan) Non-Sterile |

<b>What's Included</b>			
	<p>Easy Catcher Disposable STERILE (Case of 10)                  #O-EC51 (US), #O-EC51-EU (EU)                  #O-EC51-UK (UK), #O-EC51-JP (Japan)                  Specify Rail Type</p>		<p>Easy Catcher Disposable NON-STERILE (Case of 10)                  #O-EC52 (US), #O-EC52-EU (EU)                  #O-EC52-UK (UK), #O-EC52-JP (Japan)                  Specify Rail Type</p>
<b>What's Optional</b>	<div style="display: flex; align-items: center;">  <p>Allen Uro Collection Container, 20L with Solidifier (Case of 10)                  #O-UCC20S</p> </div>		

