

U-Step 2 Walking Stabilizer

Specifically for Neurological Conditions

Stability...

Our patented U-shaped base is ultra-stable. It surrounds you and braces you in every direction from falling. You feel as stable as you would feel holding onto another person's arm. Other walkers only partially brace you from falling.

Maneuverability...

By placing the main turning wheels towards the middle of the base, as opposed to in the rear as in other rolling walkers, the U-Step 2 turns around in 29 inches as compared to 50 inches for most other advanced walkers. This will allow you to turn in tight indoor areas, such as washrooms, kitchens, and hallways. At 23 inches wide, you can easily walk through narrow doorways.

Control...

The advanced control features play an important role in the security that the U-Step provides you. The **reversed braking system** and the **rolling resistance control** are exclusive to the U-Step 2 and put you in complete control.

Additional Unique U-Step 2 Capabilities

- Spring-Loaded Front Wheel
- Easily Transported
- Manages Curbs
- Optional Weights



Optional Feature: Laser & Sound Cueing Module

Primarily used by those with Parkinson's freezing but can be used by anyone with an irregular gait pattern. The Laser and Sound Cueing Module can help get you started, normalize your walking and increase your stride. Press the red button on the module to project a bright red laser line on the floor to guide your step. You can also turn on the sound feature to set a beat pattern for walking speed.

Medicare & Insurance Reimbursable

“When we bought a U-Step walker the results were nothing short of a miracle. Dad said he “felt safe” immediately, as if the walker were a natural extension of his body. It is lightweight, easy to operate, and with a far more dignified design than other walkers he is used to.”

— Ralph Bakery,
Los Angeles, CA

In-Step Mobility Products Corp. 1-800-558-7837

8027 Monticello Ave., Skokie, IL 60076 • Fax #: (847) 676-1202
www.ustep.com • e-mail: walkers@ustep.com • See us on Facebook and Twitter

Addressing your walking needs...

U-Step 2 Advanced features

Specifications for the U-Step 2 Walking Stabilizer

Model number	US-PC2
Medicare code (HCPC)	#E0147
Weight capacity.....	375 lbs
Height adjustment range	4'10" to 6'2"
Height customization	available
Size of Padded Seat	19" x 8"
Height from floor to seat	22"
Overall width.....	23"
Length	25"
Folded dimensions.....	approx. 42" x 23" x 10"
Turning circle	29"
Weight	20 lbs
Material of frame	Tubular Steel and Aluminum

— Patent Pending —



“ Anyone with walking difficulties should have the U-Step! Doctors and nurses have been impressed with it. My husband demonstrated it at a Parkinson’s group meeting. We are very thankful for this walker. My husband wouldn’t be able to walk without it. ”

— K. Aumann



Folds for Transport

Quality Assurance

In an effort to create the best product possible, we used only quality components to insure the U-Step 2 would hold up well over time. For example, the ergonomically designed and positioned handle bars reduce stress on the back and joints, and encourage an upright stance. The base of the unit is completely welded, not bolted, so it will not loosen over time.



- Comfort grips
- Ergonomically Positioned Handlebars
- Adjustable backrest
- Padded seat
- Height adjustable
- Laser & Sound Cueing Module (Optional)
- Removable basket
- Spring-loaded front wheel
- Place to step for going up curbs
- 4-inch non-marking casters
- Laser projected red line
- Glow-in-the-Dark tabs
- Tension control

Why the U-Step was developed

“In 1991, my mother’s neurological condition affected her walking to the point of requiring a walking aid. Her physician prescribed a standard walker, in addition to therapy and leg braces. However, the walker hindered her walking, and more importantly did not prevent her from falling. For two years, I worked with physicians, therapists and engineers to develop a better solution. Today, tens of thousands of people with neurological conditions walk more safely and easily with the U-Step” says Jonathan Miller, President of In-Step Mobility Products Corp.

Stop falling! Order your U-Step! Call 1-800-558-7837

Two week money back guarantee! • Medicare & Insurance Reimbursable